

Push-To-Set™ Adult Digital Intermittent Suction Unit (ISO/CE)

DESCRIPTION AND APPLICATION

The Push-To-Set™ Digital Intermittent Suction Unit (DPTS-ISU) exceeds clinical expectations in medical suction procedures, by offering technology breakthroughs plus the quality and simplicity you are accustomed to receiving from Ohio Medical's products. It retains the unique, patented unilogic module which provides a reliable, quiet means of cycling the regulator on and off during intermittent suctioning. It also offers dual modes to provide the option of an intermittent mode for nasogastric drainage and a continuous mode for tracheal, pharyngeal or surgical suctioning.

The PTS-ISU has an ergonomic design that contains two exclusive safety features that simplify processes, minimize clinical errors and increases product performance. The patented Push-To-Set™ mechanism automatically establishes a vacuum limit with each vacuum level setting in both continuous and intermittent modes. A unique dual-spring regulator module ensures precision in the critical care range (0-26.7 kPa/0-200 mmHg) while also providing fast adjustment (2 turns) up to full wall vacuum for emergency resuscitation procedures.

FEATURES AND BENEFITS

- **EXCLUSIVE!** Patented Push-To-Set™
Improves patient safety
May prevent inadvertent unregulated suction
Simple, one handed operation
- **EXCLUSIVE!** Quick-to-Max
Faster adjustment (2 turns from zero) to full available wall vacuum for emergency situations
- **Unilogic Timing Module**
Pneumatic, quiet, reliable operation
Returns to atmospheric pressure during "OFF" cycle
Timing cycles
Independently adjustable 3 to 30 seconds "ON" and "OFF" cycles
Factory preset at 15 seconds on/8 seconds off (± 3 seconds)
Begins in the "ON" mode
- **Flow Rates**
(Intermittent) Preset at 8 L/min (ASTM)
(Continuous) 0 to 80 L/min (without fittings at full increase)
- **Large, Easy to Read Gauge**
Color coded ranges:
Digital: Unmatched accuracy for greater safety: ± 1% of full scale
- **Safety Features**
Push-to-Set™ automatic occlusion design
Positive pressure safety relief valve
- **Modular Components**
Common components reduce service part inventory
Simplifies service procedures
Low cost replacement parts and kits
- **Durable Design and Construction**
Few mechanical parts
Rugged, shatter-resistant case
Corrosion & Lubrication Free
Service-free backplate
- **CE marked**
- **Ordering Options**
Available in various configurations including adapters for connection to source vacuum. Also available with options for collection devices and overflow safety protection.



TECHNICAL SPECIFICATIONS*

Safety Features

Push-To-Set™ Automatic occlusion device
Quick-to-Max (2 turns to full wall vacuum)

Flow Rate

Intermittent: 0-16 L/min (preset at 8 L/min per ASTM)
Continuous: 0-80 L/min (not adjustable) without fittings
at full increase

Timing

Easily adjustable Silent ON/OFF
Preset 15 seconds ON; 8 seconds OFF (±3 seconds)
Begins in the "ON" mode

Vacuum Range

0 to 26.7 kPa (0 to 200 mmHg and full wall vacuum)

Gauge Accuracy

±1% of full scale (2.0 mmHg)

Dimensions

16.5 x 7.1 x 12.2 cm (6.5" H X 2.8" W X 4.8" D)

Weight: 0.57 Kg (1 lb 4oz) less fittings and adapters

F.O.B. Point: Gurnee, IL

* Specifications are nominal, subject to change without notice. Please see "Instructions for use" manual for the most current regulated vacuum ranges.

Applied Standards

ISO 10079-3 and ASTM F960

SHIPPING INFORMATION

Shipping Weight: 0.77 Kg (1 lb 11 oz)
less fittings and adapters

Package Size: 16.5 x 26.67 x 20.32 cm
(6.5" H x 10.5" W x 8" D)

Packaging: One per carton

Terms: net 30 days



Parts Configuration

87XX - XXXX - 9XX

Push-To-Set™ Series

Adapter (vacuum source connection)

No Adapter with 1/8" NPT Female.....	00
BOC Type (British Standard).....	31
DIN Type (German Standard).....	32
AFNOR Type (French Standard).....	33
AGA Type (Nordic Standard).....	34

Fittings (patient connection)

No Fitting 1/8" NPT Female.....	00
Tubbing Nipple.....	01
Locking Gland.....	04
Locking Gland & Overflow Safety Trap.....	05
Locking Gland, Overflow Safety Trap, 1300 ml Plastic Bottle with Wall Slide, Tubing.....	06

Push-To-Set™ Digital Intermittent Suction Unit

1353	Off-White*	1753	Red
1453	Yellow	1853	Lt. Blue
1553	Blue	1953	Purple
1653	Pink		

*Standard Color

© 2017 Ohio Medical LLC. All rights reserved.

This document contains information that is proprietary and confidential to Ohio Medical LLC. Use of this information is under license from Ohio Medical LLC. Any use other than that authorized by Ohio Medical LLC is prohibited.



Ohio Medical and the Ohio Medical logo are registered trademarks and
Push-To-Set is a trademark of Ohio Medical LLC.
U.S. Patent #6,623,463 (PushTo-Set™ Vacuum System)

Ohio Medical
1111 Lakeside Drive - Gurnee, IL 60031 USA
Phone: +1 847 855 0500 - Fax: +1 847 855 6218
www.ohiomedical.com
Form No. 255451 (Rev. 5, 01/2017)

Ohio Medical Authorized Representative



OxygenCare Ltd.
2 Holfeld Business Park
Kilmacanogue Co Wicklow Ireland
Tel: +353 1 2769700 - Fax: +353 1 2764970