

Intermittent Suction Unit (ANSI) Adult

DESCRIPTION AND APPLICATION

The Ohio Medical® Intermittent Suction Unit (ISU) is a combination continuous and intermittent regulator for use with hospital pipeline vacuum systems. It provides the versatility of an intermittent mode for nasogastric drainage and a continuous mode for tracheal, pharyngeal, or surgical suctioning. In the intermittent mode, the unit uses a pneumatic logic device to alternately turn the suction “ON” and “OFF.” This unique Unilogic Module provides quiet and reliable controlled intermittent suction. The unit is equipped with independent controls to adjust the duration of the “ON” and “OFF” cycles, the level of vacuum, and the rate of flow.

The regulator also provides continuous vacuum, adjustable over the full range of available vacuum, to meet all clinical continuous suction requirements.

The Intermittent Suction Unit is classified as a high vacuum/high flow regulator.



Adult Intermittent Suction Unit

FEATURES AND BENEFITS

- **Unilogic Module**
 - Pneumatic
 - Quiet, reliable operation
 - Returns to atmospheric pressure during “OFF” cycle
- **Fully Adjustable**
 - Timing cycles
 - 3 to 30 seconds
 - “ON” and “OFF” independently adjustable
 - Factory preset at 15 seconds on/8 seconds off (begins in the off mode)
- **Flow rates**
 - (Intermittent) 0 to 13 LPM, preset at 8 LPM (ASTM)
 - (Continuous) 0 to 80 LPM (not adjustable)
- **Large, Easy to Read Gauge**
 - Color coded ranges
 - 5 mmHg increments
 - Glow-in-the-dark gauge increments and needle
 - Accuracy $\pm 5\%$ full scale
- **Modular Components**
 - Simplifies service procedures
 - Low cost replacement parts and kits
- **Durable Design and Construction**
 - Few mechanical parts
 - Rugged, shatter-resistant case
- **Ordering Options**
 - Available in various configurations to include adapters for connection to source vacuum. Also available with options for collection devices and overflow safety protection.

SPECIFICATIONS*

Flow Rate

Intermittent : 0-13 LPM (preset at 8 LPM per ASTM)
 Continuous: 0-80 LPM (not adjustable) without fittings at full increase

Timing

3 to 30 seconds
 Independent on and off
 Preset 15 seconds on/8 seconds off
 Begins in the "OFF" mode

Vacuum Range

0 to full range of available vacuum

Gauge Range

0-200 mmHg and full vacuum

Gauge Accuracy

±10 mmHg or ±5% Full Scale Deflection

Dimensions

6.6" H x 3.5" W x 4.8" D
 (16.8 x 8.9 x 12.2 cm)

Weight

1 lb 11oz (0.781 Kg)

F.O.B. Point

Gurnee, IL

Applied Standard

ISO 10079-3

* Specifications are nominal, subject to change without notice.

SHIPPING INFORMATION

Shipping Weight

2 lbs, 7 oz. (1.1 kg)

Package Size

8" H x 3.5" W x 4.7" D
 (20.3 x 8.9 x 12 cm)


Packaging

One per carton

Terms

Net 30 days

Parts Configuration



© 2018 Ohio Medical, LLC.

This document contains information that is proprietary and confidential to Ohio Medical, LLC. Use of this information is under license from Ohio Medical, LLC. Any use other than that authorized by Ohio Medical, LLC is prohibited.

Ohio Medical and the Ohio Medical logo are registered trademarks of Ohio Medical, LLC.

Puritan Bennett is a registered trademark of Covidien LP.
 Chemetron is a registered trademark of Allied Healthcare Products, Inc.



Ohio Medical LLC

1111 Lakeside Drive, Gurnee, IL 60031 USA
 Phone: 1.866.549.6446
www.ohiomedical.com
 550043-ISU (Rev. 3) 04/2018