

Thoracic Vacuum Regulator (ISO)

DESCRIPTION AND APPLICATION

The Ohio Medical® Thoracic Vacuum Regulator brings precision and safety to your pleural, mediastinal and other low suction applications. It is designed to maintain accurate and stable negative pressures at low levels of suction while maximizing flow rates. It compensates for positive pressure in the patient line (e.g. patient coughing) with a positive pressure relief valve which immediately opens to vent the pressure. The Thoracic vacuum regulator also provides high flow response to air leaks in the lung. It adds convenience, flexibility, and reliability to disposable or reusable chest drainage systems, and can add cost savings when used in conjunction with disposables. The regulator must be used with an underwater seal and collection bottle (disposable or reusable). It is a silent alternative to noisy bubbling water manometers.

The Thoracic Vacuum Regulator is available with probes for virtually all pipeline systems and is offered with a standard tubing nipple.

FEATURES AND BENEFITS

- **Automatically compensates for vacuum fluctuations** at the catheter (e.g., air leaks) or wall outlet.
- **Large, precision gauge** with color coded needle to indicate low vacuum.
- **Fine-tuning suction control knob** moves slowly and precisely, taking 5 full turns to span the entire gauge range.
- **Negative pressure may be increased or decreased quickly** and easily without the use of water, funnels, or syringes.
- **Positive pressure relief valve** vents positive pleural pressure or accidental positive pressure away from the patient and system.
- **Separate toggle ON/OFF** switch retains previous settings.
- **Extended range** provides higher negative pressures required for multiple chest tubes and other applications.
- **Pressure limiting valve** prevents negative pressure from exceeding 55 cm H₂O (5.4 kPa).
- **Delivers high flow rates** at least 40 Lpm at full increase dependent upon supply vacuum and setup.
- **Silent vacuum control** eliminates bubbling noise of water manometers.
- **Full Three Year Parts and Labor Warranty**
- **CE marked**
- **Ordering Options**
Available in various configurations to include Probes for connection to source vacuum. Also available with options for collection devices and overflow safety protection.



Thoracic Vacuum Regulator

- 0-65 cm H₂O
- 0-6 kPa

Specifications*

Flow Rate

0-40 Lpm at full increase

Gauge Accuracy

± 3 cm H₂O

Dimensions

18.29 x 8.9 x 10.92 cm
(7.2" x 3.5" x 4.3")

Weight

0.539 Kg (19 oz) (less fittings)

Warranty

Three Years

Applied Standard

ISO 10079-3

* Specifications are nominal, subject to change without notice.

Purchasing Specifications

Thoracic Vacuum Regulator with color coded gauge marked in 1 cm H₂O increments. Regulator provides continuous, regulated suction from a piped vacuum system. Regulator can be used in conjunction with disposable or reusable chest drainage systems. Adjustable, pressure limiting valve prevents pressures from exceeding set level. Positive pressure relief valve vents positive pleural pressures or accidental pressures in excess of 10 cm H₂O (1.0 kPa) away from patient system. Separate on/off switch conserves vacuum, allows retention of previous setting.

Shipping Information

Shipping Weight

0.71 Kg (1 lb 9 oz) less fittings and adapters

Package Size

26.16 H x 19.8 W x 17.27 D cm
(10.3" H x 7.8" W x 6.8" D)

Packaging

One per carton

Parts Configuration

67XX - XXXX - 9XX

Adapter/Probe (vacuum source connection)

No Adapter with 1/8" NPT Female.....	00
British Standard (BOC).....	31
German Standard (DIN).....	32
French Standard (AFNOR).....	33
Nordic Standard (AGA).....	34
Australian (SIS).....	35
Japanese.....	14

Fittings (patient connection)

No Adapter with 1/8" NPT Female.....	00
Tubing Nipple.....	01
Locking Gland and Overflow Safety Trap.....	05

Thoracic Vacuum Regulator

English.....	1276
German.....	1323
Dutch.....	1385
Spanish.....	1278
French.....	1277
Greek.....	1384
Italian.....	1380
Japanese.....	1382
Portuguese.....	1386
Russian.....	1383
Swedish.....	1381

© 2017 Ohio Medical LLC.

This document contains information that is proprietary and confidential to Ohio Medical LLC. Use of this information is under license from Ohio Medical LLC. Any use other than that authorized by Ohio Medical LLC is prohibited.

Ohio Medical and the Ohio Medical logo are registered trademarks of Ohio Medical LLC.



Ohio Medical
1111 Lakeside Drive - Gurnee, IL 60031 USA
Phone: +1 847 855 0500 - Fax: +1 847 855 6218
www.ohiomedical.com
SOT122 (Rev. 7, 01/2017)



Ohio Medical Authorized Representative



OxygenCare Ltd.
2 Holfeld Business Park
Kilmacanogue Co Wicklow Ireland
Tel: +353 1 2769700 - Fax: +353 1 2764970