


# Wall/Console Slide Bracket & Blank Wall Assembly



## PART NUMBERS

- 261690 Wall Slide
- 261691 Console Slide
- 261394 Wall Cover Blank

## MATERIALS OF CONSTRUCTION

- Slide: Aluminum
- Screws: Plated Steel
- Cover Plate: Cycloyol®

Cycloyol® is a registered trademark of GE Plastics.

Ohio Medical Corporation  
1111 Lakeside Drive  
Gurnee, IL 60031-4099 USA  
800-448-0770 • www.ohiomedical.com  
Ph: 847-855-0500 • Fax: 847-855-6300

255093 (Rev.3) 02/06

