

Hose Assemblies and Hose Retractor



SPECIFICATIONS

The Ohio Medical Healthcair® Series hose assemblies are manufactured with FDA approved medical grade anti-static PVC tubing. The tubing is reinforced with polyester braiding. The maximum working pressure is 200 psi @ 70°F. The fittings have either the coupling or the grips color coded per CGA C-9 standards. The fittings are permanently crimped on each end to the tubing with ferrules.

The hose assemblies are normally made with a DISS (Diameter Index Safety System) female connection using a handtight nipple with a wrench tight nut on one end (to be connected to an outlet). The other end is a male DISS (Diameter Index Safety System) or female quick connect coupler (Ohmeda or Chemetron) with check unit. The standard length is approximately 5 or 6 feet. The internal diameters on Ohio Medical Healthcair® series hose assemblies are 1/4" for pressurized gases and 5/16" for vacuum & WAGD hoses.

The Ohio Medical Healthcair® Series hose retractor can be mounted directly on the ceiling outlet or to the ceiling. The retractor wire is manufactured out of 48" stainless steel heavy-duty cable.

PART NUMBERS

DISS Female to Male DISS Male - 5' Hose Length

263400-1-05	Oxygen Hose Assembly
263400A-2-05	Vacuum Hose Assembly
263400-3-05	Nitrous Oxide Hose Assembly
263400-4-05	Air Hose Assembly
263400-5-05	Nitrogen Hose Assembly
263400A-6-05	WAGD Hose Assembly
263400-7-05	Carbon Dioxide Hose Assembly

DISS Female to Ohmeda Compatible Quick Connect - 5' Hose Length

263401-1-05	Oxygen Hose Assembly
263401A-2-05	Vacuum Hose Assembly
263401-3-05	Nitrous Oxide Hose Assembly
263401-4-05	Air Hose Assembly
263401A-6-05	WAGD Hose Assembly

DISS Female to Chemetron Compatible Quick Connect - 5' Hose Length

263402-1-05	Oxygen Hose Assembly
263402A-2-05	Vacuum Hose Assembly
263402-3-05	Nitrous Oxide Hose Assembly
263402-4-05	Air Hose Assembly
263402A-6-05	WAGD Hose Assembly

Retractor

261746	Retractor
--------	-----------

TRUSTED BRANDS OF OHIO MEDICAL



Ohio Medical Corporation - 1111 Lakeside Drive - Gurnee, IL 60031-4099 USA
 TollFree: 800-448-0770 - Phone: 847-855-0500 - Fax: 847-855-6300
www.ohiomedical.com