



Ohio
Medical Corporation®

always there | for life





always there | for life



Ohio Medical Corporation

has been the leading supplier of Healthcair® medical air and vacuum pumping systems, Ohio Medical Corporation® and Amvex medical suction regulators, Ohio Medical™ and Amvex air and oxygen flow meters and Aeros® portable suction equipment since 1958.

Ohio Medical Corporation is an ISO 13485 registered and FDA regulated medical device manufacturing facility. All Ohio Medical medgas products are produced to conform to NFPA 99 standards.

With a leading research and development team, dedicated service department, and knowledgeable customer support teams, we will continue to provide the medical industry with leading technology and quality products as we have for the past five decades.

Ohio Medical is the only name you need to know for your medical suction and oxygen therapy and medgas equipment needs.







your patients|your partner



Ohio Medical

has maintained a presence in the medical industry for well over 45 years under the Ohio Medical, Amvex, and Aeros-Instruments & Squire-Cogswell names. We are a leading supplier of medical pumping systems, with a sound reputation in hospital environments around the world.

In the 1950 the first Ohio Medical Vacuum regulator was introduced under the name of Ohio Chemical, followed by the Intermittent Vacuum regulator in 1951 and the HH Model in 1953. These regulators were manufactured by Squire-Cogswell (now Ohio Medical) and private labeled for Ohio Medical, Ohmeda and successive companies since their introduction over 45 years ago.

In July of 2005, Ohio Medical Corporation was formed, combining the Ohmeda® Suction & Oxygen Therapy products, Aeros portable suction equipment, the Healthcair and Selectorr pumping systems and medgas pipeline products. Finally in November 2007 Ohio Medical acquired Amvex Corporation a leading supplier of medical gas devices. This union created a synergistic blend of products and technology with strong customer and patient safety focus.







the priority | of quality



Quality is essential


At Ohio Medical, we integrate quality into all aspects of our business, from product development to customer support. Our high level of quality comes from our commitment to two key issues:

- Meeting or surpassing the regulatory standards in the countries where we make and sell our products.
- Maintaining a staff of talented employees whose commitment to quality supports Ohio Medical's focus.

Ohio Medical Corporation is an ISO 13485 registered and FDA regulated medical device manufacturing facility.

A major part of Ohio Medical's quality efforts entail certification according to domestic and global standards such as UL, CSA, NFPA, NEC, ISO 13485, CMDCAS, ETL and CE through the Medical Device Directive (MDD).







clinical | solutions



Suction Regulators

Ohio Medical offers a wide selection of suction regulators for the industry. Our regulators are specially designed for a variety of clinical applications including airway, nasogastric and surgical suctioning as well as chest drainage. Ohio Medical has been providing suction products to health-care facilities worldwide for over 45 years.

Oxygen Therapy/ Respiratory Care

Ohio Medical's Pressure Compensated Flowmeters are designed to meet strict standards of durability and precision. Several models for delivery of oxygen and air are available to meet your needs.

Portable & Venturi Suction

A full line of self-standing and table top suction machines is available for procedures that require intermittent or constant suction. Both AC and battery powered models are offered.

The Ohio Medical Venturi gas-powered portable vacuum unit provides a convenient suction source on infant radiant warmers, during transport, in special procedure areas or on crash carts. When powered by oxygen, the Du-O-Vac Plus also provides two O₂ DISS outlets for oxygen therapy applications.

Gas Fittings & Connectors

Customer convenience is assured by our complete line of gas fittings and connectors. With models to meet every need, they are reliable, easy to use, safe and familiar.







system | solutions



Vacuum Systems

Rotary Vane, Rotary Claw, Liquid Ring

Our vacuum systems comply with the latest NFPA 99 standards. We offer compact modular tank and stack mount systems, in standard and custom configurations. Our choices include rotary vane, liquid ring and rotary claw.

Air Systems

Oil-Less Scroll, Oil-Less Reciprocating

Oil-less breathing air systems are offered in both piston and scroll versions. All models are NFPA 99 compliant with all needed equipment, such as CO/dewpoint monitor, dryers, filters, etc.

Medgas Pipeline

We offer a full line of NFPA 99 compliant medgas pipeline products including outlets in a variety of connection styles, area and master alarms, manifolds, valve boxes, and many other products to meet all of your needs.

Vacuum Pumps

Let us help with all of your vacuum pump replacements. Whether your system is equipped with rotary vane or liquid ring type pumps, we can help. Also, contact us about upsizing your current system.

Air Compressors

Our lines of oil-less air compressors are available in both reciprocating and scroll configurations. Standard compressor modules include motors, inlet filters, belts, belt guards, and aftercoolers.

Accessories

We inventory a vast array of replacement parts for many makes of air and vacuum systems. Call us for oil, filters, gauges, relief valves, dryers, tanks, etc.





1111 Lakeside Drive, Gurnee, IL 60031-4099 USA | 1.800.448-0770 | www.ohiomedical.com