


Adapters, Couplers, & Hose Assemblies



Bulk Hoses

- Available in custom lengths and configurations
- Available in a variety of gases
- Available in USA and ISO colors



Gas	Conductive	Non-Conductive
OXYGEN	HS-01U-OC4	HS-01U-ON4
MEDICAL AIR	HS-01U-AC4	HS-01U-AN4
VACUUM 1/4"	HS-01U-VC4	HS-01U-VN4
VACUUM 5/16"	HS-01U-VC5	N/A
N ₂ O	HS-01U-2C4	HS-01U-2N4
CO ₂	HS-01U-CC4	HS-01U-CN4
WAGD 1/4"	HS-01U-WC4	N/A
WAGD 5/16"	HS-01U-WC5	HS-01U-WN5
NITROGEN	HS-01U-NC4	HS-01U-NN4
INSTRUMENT AIR	HS-01U-IC4	N/A
GENERAL PURPOSE	N/A	HS-01U-GN4

For ISO options, please replace "U" with "I".

Bulk Hose Accessories



Heavy Duty Hose Retractor
HS-RET-HD



Hose Clamp for 1/4" Hose
HS-CLAMP-500




Ferrules for 1/4" Hose

HS-FER-562 3/4" Long
HS-FER-562-1 1" Long

Ferrules for 5/16" Hose

HS-FER-VAC 3/4" Long
HS-FER-VAC-1 1" Long



Hose Drop Assembly Part Numbers

DISS Female to Male DISS Male - 5' Hose Length

HS-05UO-DFDCC4	Oxygen
HS-05UV-DFDCC5	Vacuum
HS-05U2-DFDCC4	Nitrous Oxide
HS-05UA-DFDCC4	Medical Air
HS-05UN-DFDCC4	Nitrogen
HS-05UW-DFDCC5	WAGD
HS-05UC-DFDCC4	Carbon Dioxide

Chemetron Quick Connects - 5' Hose Length

HS-05UO-CFDFC4	Oxygen
HS-05UV-CFDFC5	Vacuum
HS-05U2-CFDFC4	Nitrous Oxide
HS-05UA-CFDFC4	Medical Air
HS-05UW-CFDFC5	WAGD

Diamond/Ohmeda Quick Connects - 5' Hose Length


HS-05UO-DFOFC4	Oxygen
HS-05UV-DFOFC5	Vacuum
HS-05U2-DFOFC4	Nitrous Oxide
HS-05UA-DFOFC4	Medical Air
HS-05UW-DFOFC5	WAGD











Hose Drop



Hose Assemblies

- Available in custom lengths and configurations
- Available in a variety of gases
- Available in USA and ISO colors
- Hose assemblies also available with international fittings (NIST, British, etc.)



Please contact us for information on our MR Conditional options	 OHMEDA ADAPTER	 CHEMETRON ADAPTER	 DISS HANDTIGHT	 DISS HEX NUT	 DISS MALE WITH DEMAND CHECK
 OHMEDA COUPLER	HS-LLCG-OFOMC4	HS-LLCG-OFCMC4	HS-LLCG-DHOF C4	HS-LLCG-DFOFC4	HS-LLCG-OFDCC4
 CHEMETRON COUPLER	HS-LLCG-CFOMC4	HS-LLCG-CFCMC4	HS-LLCG-CFDHC4	HS-LLCG-CFDFC4	HS-LLCG-CFDCC4
 DISS HANDTIGHT	HS-LLCG-DHOMC4	HS-LLCG-DHCMC4	HS-LLCG-DHDHC4	HS-LLCG-DFDHC4	HS-LLCG-DHDCC4
 DISS HEX NUT	HS-LLCG-DFOMC4	HS-LLCG-DFCMC4	HS-LLCG-DFDHC4	HS-LLCG-DFDFC4	HS-LLCG-DFDCC4
 DISS MALE WITH DEMAND CHECK	HS-LLCG-DCOMC4	HS-LLCG-CMDCC4	HS-LLCG-DHDCC4	HS-LLCG-DFDCC4	HS-LLCG-DCDCC4

Part Configurator

HS-LLCG-XXYYZZ(U)

Length (ft)					
Color Code	Gas Type	Hose Type	Options		
USA U	Nitrous Oxide 2	Conductive 1/4" C4	MRI M		
ISO I	Medical Air A	Conductive 5/16" C5	Brass Fitting B		
Blue B	Carbon Dioxide C	DEHP Free 1/4" D4	DIN Rail Clamp C		
Black Non-Conductive K	Heliox H	General Purpose 1/4" G4	DIN Rail Slide D		
Pin-Stripe (ISO Air) P	Instrument Air I	Non-Conductive 1/4" N4			
Stripe (ISO Air) S	O ₂ He J	Non-Conductive 5/16" N5			
White W	Nitrogen N				
	Oxygen O				
	Medical Vacuum V				
	WAGD W				














For correct hose matrix, follow these rules when selecting inlet and outlet connections:

- Female fittings first
- Numbers before letters
- In alphabetical order

Inlet (XX) & Outlet (YY) Connections									
Blank XX	3/8" FNPT F6	AGA Female ZF	French Male Angled FE	Ohmeda Male OM					
Female Fittings	DISS Handtight with Elbow HE	AGA Female with Elbow ZE	French Male FM	Puritan Male with Elbow PE					
1/4" FNPT Brass Swivel 4S	Japanese Female JF	Male Fittings	German Male no Check Angled GE	Puritan Male PM					
Australian Handtight AH	MedStar Female MF	M16 Air Adapter 16	German Male Straight GM	Carburos Male with Elbow RE					
British Female BF	NIST Nut with Elbow NE	1/8" MNPT with Elbow 2E	Japanese Male JM	Carburos Male RM					
Chemetron Female CF	NIST Nut Female NF	Australian Male AM	1/8" MNPT M2	Schrader Male SM					
DISS Hex Nut with Elbow DA	Ohmeda Female OF	British Male BM	1/4" MNPT M4	Schrader Swivel Adapter SS					
DISS Hex Nut DF	Pole Clamp F2 1.5" PC	Chemetron Male with Elbow CE	3/8" MNPT M6	1/4" Hose Barb Tee T4					
DISS Handtight DH	Puritan Female PF	Chemetron Male CM	MedStar Male MM	Oxequip Male with Elbow XE					
1/8" FNPT F2	Pole Clamp F2 2.5" Long PL	DISS Male with Demand Check DC	NIST Male Angled NA	Oxequip Male XM					
1/4" FNPT F4	Carburos Female RF	DISS Male no Demand Check with Elbow DE	NIST Male NM	AGA Male ZM					
	Schrader Female SF	DISS Male no Demand Check DM	Ohmeda Male with Elbow OE						











NOTE: Not all configurations are possible. Other options are available, please contact your Sales Representative for more information.

Ohmeda® Adapters

Please contact us for information on our MR Conditional options						
	OXYGEN	MEDICAL AIR	VACUUM	N₂O	WAGD	CO₂
 1/8" NPT MALE	AD-OHUO-M2	AD-OHUA-M2	AD-OHUV-M2	AD-OHU2-M2	AD-OHUW-M2	AD-OHUC-M2
 1/4" HOSE BARB 5/16" HOSE BARB	AD-OHUO-HB4	AD-OHUA-HB4	AD-OHUV-HB4 AD-OHUV-HB5	AD-OHU2-HB4	AD-OHUW-HB4 AD-OHUW-HB5	AD-OHUC-HB4
 DISS MALE	AD-OHUO-DM	AD-OHUA-DM	AD-OHUV-DM	AD-OHU2-DM	AD-OHUW-DM	AD-OHUC-DM
 1/4" NPT MALE	AD-OHUO-M4	AD-OHUA-M4	AD-OHUV-M4	AD-OHU2-M4	AD-OHUW-M4	AD-OHUC-M4
 1/8" NPT FEMALE	AD-OHUO-F2	AD-OHUA-F2	AD-OHUV-F2	AD-OHU2-F2	AD-OHUW-F2	AD-OHUC-F2
 DISS HANDTIGHT	AD-OHUO-DH	AD-OHUA-DH	AD-OHUV-DH	AD-OHU2-DH	AD-OHUW-DH	AD-OHUC-DH
 DISS NUT & NIPPLE	AD-OHUO-DN	AD-OHUA-DN	AD-OHUV-DN	AD-OHU2-DN	AD-OHUW-DN	AD-OHUC-DN














For ISO colors, replace the "U" with "I" for Oxygen, Medical Air and Vacuum only.

Ohmeda® Couplers

<p>Please contact us for information on our MR Conditional options</p>						
	OXYGEN	MEDICAL AIR	VACUUM	N₂O	WAGD	CO₂
 1/8" NPT MALE	CO-OHUO-M2	CO-OHUA-M2	CO-OHUV-M2	CO-OHU2-M2	CO-OHUW-M2	CO-OHUC-M2
 1/4" HOSE BARB 5/16" HOSE BARB	CO-OHUO-HB4	CO-OHUA-HB4	CO-OHUV-HB4 CO-OHUV-HB5	CO-OHU2-HB4	CO-OHUW-HB4 CO-OHUW-HB5	CO-OHUC-HB4
 DISS MALE	CO-OHUO-DM	CO-OHUA-DM	CO-OHUV-DM	CO-OHU2-DM	CO-OHUW-DM	CO-OHUC-DM
 1/4" NPT MALE	CO-OHUO-M4	CO-OHUA-M4	CO-OHUV-M4	CO-OHU2-M4	CO-OHUW-M4	CO-OHUC-M4
 1/8" NPT FEMALE	CO-OHUO-F2	CO-OHUA-F2	CO-OHUV-F2	CO-OHU2-F2	CO-OHUW-F2	CO-OHUC-F2
 DISS HANDTIGHT	CO-OHUO-DH	CO-OHUA-DH	CO-OHUV-DH	CO-OHU2-DH	CO-OHUW-DH	CO-OHUC-DH
 DISS NUT & NIPPLE	CO-OHUO-DN	CO-OHUA-DN	CO-OHUV-DN	CO-OHU2-DN	CO-OHUW-DN	CO-OHUC-DN














For ISO colors, replace the "U" with "I" for Oxygen, Medical Air and Vacuum only.

Chemetron® Adapters

<p>Please contact us for information on our MR Conditional options</p>						
	OXYGEN	MEDICAL AIR	VACUUM	N₂O	WAGD	CO₂
 1/8" NPT MALE	AD-CHUO-M2	AD-CHUA-M2	AD-CHUV-M2	AD-CHU2-M2	AD-CHUW-M2	AD-CHUC-M2
 1/4" HOSE BARB 5/16" HOSE BARB	AD-CHUO-HB4	AD-CHUA-HB4	AD-CHUV-HB4 AD-CHUV-HB5	AD-CHU2-HB4	AD-CHUW-HB4 AD-CHUW-HB5	AD-CHUC-HB4
 DISS MALE	AD-CHUO-DM	AD-CHUA-DM	AD-CHUV-DM	AD-CHU2-DM	AD-CHUW-DM	AD-CHUC-DM
 1/4" NPT MALE	AD-CHUO-M4	AD-CHUA-M4	AD-CHUV-M4	AD-CHU2-M4	AD-CHUW-M4	AD-CHUC-M4
 1/8" NPT FEMALE	AD-CHUO-F2	AD-CHUA-F2	AD-CHUV-F2	AD-CHU2-F2	AD-CHUW-F2	AD-CHUC-F2
 DISS HANDTIGHT	AD-CHUO-DH	AD-CHUA-DH	AD-CHUV-DH	AD-CHU2-DH	AD-CHUW-DH	AD-CHUC-DH
 DISS NUT & NIPPLE	AD-CHUO-DN	AD-CHUA-DN	AD-CHUV-DN	AD-CHU2-DN	AD-CHUW-DN	AD-CHUC-DN









For ISO colors, replace the "U" with "I" for Oxygen, Medical Air and Vacuum only.

Chemetron® Couplers









      <p>Please contact us for information on our MR Conditional options</p>						
	OXYGEN	MEDICAL AIR	VACUUM	N₂O	WAGD	CO₂
 <p>1/8" NPT MALE</p>	CO-CHUO-M2	CO-CHUA-M2	CO-CHUV-M2	CO-CHU2-M2	CO-CHUW-M2	CO-CHUC-M2
 <p>1/4" HOSE BARB 5/16" HOSE BARB</p>	CO-CHUO-HB4	CO-CHUA-HB4	CO-CHUV-HB4 CO-CHUV-HB5	CO-CHU2-HB4	CO-CHUW-HB4 CO-CHUW-HB5	CO-CHUC-HB4
 <p>DISS MALE</p>	CO-CHUO-DM	CO-CHUA-DM	CO-CHUV-DM	CO-CHU2-DM	CO-CHUW-DM	CO-CHUC-DM
 <p>1/4" NPT MALE</p>	CO-CHUO-M4	CO-CHUA-M4	CO-CHUV-M4	CO-CHU2-M4	CO-CHUW-M4	CO-CHUC-M4
 <p>1/8" NPT FEMALE</p>	CO-CHUO-F2	CO-CHUA-F2	CO-CHUV-F2	CO-CHU2-F2	CO-CHUW-F2	CO-CHUC-F2
 <p>DISS HANDTIGHT</p>	CO-CHUO-DH	CO-CHUA-DH	CO-CHUV-DH	CO-CHU2-DH	CO-CHUW-DH	CO-CHUC-DH
 <p>DISS NUT & NIPPLE</p>	CO-CHUO-DN	CO-CHUA-DN	CO-CHUV-DN	CO-CHU2-DN	CO-CHUW-DN	CO-CHUC-DN

For ISO colors, replace the "U" with "I" for Oxygen, Medical Air and Vacuum only.

Puritan Bennett® Adapters and Couplers





Please contact us for information on our MR Conditional options	OXYGEN	MEDICAL AIR	VACUUM	N ₂ O	WAGD
  1/8" NPT MALE	AD-PBUO-M2 CO-PBUO-M2	AD-PBUA-M2 CO-PBUA-M2	AD-PBUV-M2 CO-PBUV-M2	AD-PBU2-M2 CO-PBU2-M2	AD-PBUW-M2 CO-PBUW-M2
  1/4" HOSE BARB AND 5/16" HOSE BARB	AD-PBUO-HB4 CO-PBUO-HB4	AD-PBUA-HB4 CO-PBUA-HB4	AD-PBUV-HB4 AD-PBUV-HB5 CO-PBUV-HB4 CO-PBUV-HB5	AD-PBU2-HB4 CO-PBU2-HB4	AD-PBUW-HB4 AD-PBUW-HB5 CO-PBUW-HB4 CO-PBUW-HB5
  1/4" NPT MALE	AD-PBUO-M4 CO-PBUO-M4	AD-PBUA-M4 CO-PBUA-M4	AD-PBUV-M4 CO-PBUV-M4	AD-PBU2-M4 CO-PBU2-M4	AD-PBUW-M4 CO-PBUW-M4
  1/8" NPT FEMALE	AD-PBUO-F2 CO-PBUO-F2	AD-PBUA-F2 CO-PBUA-F2	AD-PBUV-F2 CO-PBUV-F2	AD-PBU2-F2 CO-PBU2-F2	AD-PBUW-F2 CO-PBUW-F2

Schrader® Adapters and Couplers

Please contact us for information on our MR Conditional options	OXYGEN	MEDICAL AIR	VACUUM	N ₂ O	NITROGEN
  1/8" NPT MALE	AD-SHUO-M2 CO-SHUO-M2	AD-SHUA-M2 CO-SHUA-M2	AD-SHUV-M2 CO-SHUV-M2	AD-SHU2-M2 CO-SHU2-M2	AD-SHUN-M2 CO-SHUN-M2
  1/4" HOSE BARB AND 5/16" HOSE BARB	AD-SHUO-HB4 CO-SHUO-HB4	AD-SHUA-HB4 CO-SHUA-HB4	AD-SHUV-HB4 AD-SHUV-HB5 CO-SHUV-HB4 CO-SHUV-HB5	AD-SHU2-HB4 CO-SHU2-HB4	AD-SHUN-HB4 CO-SHUN-HB4
  1/4" NPT MALE	AD-SHUO-M4 CO-SHUO-M4	AD-SHUA-M4 CO-SHUA-M4	AD-SHUV-M4 CO-SHUV-M4	AD-SHU2-M4 CO-SHU2-M4	AD-SHUN-M4 CO-SHUN-M4
  1/8" NPT FEMALE	AD-SHUO-F2 CO-SHUO-F2	AD-SHUA-F2 CO-SHUA-F2	AD-SHUV-F2 CO-SHUV-F2	AD-SHU2-F2 CO-SHU2-F2	AD-SHUN-F2 CO-SHUN-F2









For ISO colors, replace the "U" with "I". For swivel options on the adapter, add "S" to the end of the part number.

Oxequip® and MedStar® Adapters

Please contact us for information on our MR Conditional options	OXYGEN	MEDICAL AIR	VACUUM	N ₂ O
 1/8" NPT MALE OXEQUIP	AD-XQUO-M2	AD-XQUA-M2	AD-XQUV-M2	AD-XQU2-M2
 1/4" HOSE BARB OXEQUIP 5/16" HOSE BARB OXEQUIP	AD-XQUO-HB4	AD-XQUA-HB4	AD-XQUV-HB4 AD-XQUV-HB5	AD-XQU2-HB4
 1/8" NPT MALE MEDSTAR	AD-MSUO-M2	AD-MSUA-M2	AD-MSUV-M2	AD-MSU2-M2
 1/4" HOSE BARB MEDSTAR 5/16" HOSE BARB MEDSTAR	AD-MSUO-HB4	AD-MSUA-HB4	AD-MSUV-HB4 AD-MSUV-HB5	AD-MSU2-HB4

For ISO colors, replace the "U" with "I". Please contact us for additional options on our Medstar couplers.

International Fittings

 BS - British	 DIN - German	 AGA - Scandinavian	 AFNOR - French
 SIS - Australian	 UNIFOR - Italian	 Japanese - Korean	 NIST - European

For gases other than Oxygen please replace the "O" in the item number with the following: "A": Medical Air, "2": Nitrous Oxide, "V": Vacuum.
For fittings options, see chart below.

1/8" MNPT	1/4" MNPT	1/4" Hose Barb	1/8" FNPT	1/4" FNPT	5/16" Hose Barb
M2	M4	HB4	F2	F4	HB5

DISS Fittings

Please contact us for information on our MR Conditional options.		OXYGEN	MEDICAL AIR	VACUUM	N ₂ O	NITROGEN	WAGD	CO ₂	INSTRUMENT AIR
1/8" NPT MALE		BA-O-M2	BA-A-M2	BA-V-M2	BA-2-M2	BA-N-M2	BA-W-M2	BA-C-M2	BA-IA-M2
		DA-15O-M2-DN	DA-15A-M2-DN	DA-15V-M2-DN	DA-152-M2-DN	DA-15N-M2-DN	DA-15W-M2-DN	DA-15C-M2-DN	DA-15I-M2-DN
		DA-15UO-M2-DH	DA-15UA-M2-DH	DA-15UV-M2-DH	DA-15U2-M2-DH	DA-15UN-M2-DH	DA-15UW-M2-DH	DA-15UC-M2-DH	DA-15UI-M2-DH
1/4" & 5/16" HOSE BARB		BA-O-HB4	BA-A-HB4	BA-V-HB4 BA-V-HB5	BA-2-HB4	BA-N-HB4	BA-W-HB4 BA-W-HB5	BA-C-HB4	BA-IA-HB4
		DA-HB4-O-DN	DA-HB4-A-DN	DA-HB4-V-DN DA-HB5-V-DN	DA-HB4-2-DN	DA-HB4-N-DN	DA-HB4-W-DN DA-HB5-W-DN	DA-HB4-C-DN	DA-HB4-I-DN
		DA-HB4-UO-DH	DA-HB4-UA-DH	DA-HB4-UV-DH DA-HB5-UV-DH	DA-HB4-U2-DH	DA-HB4-UN-DH	DA-HB4-UW-DH DA-HB5-UW-DH	DA-HB4-UC-DH	DA-HB4-UI-DH
1/4" NPT MALE		BA-O-M4	BA-A-M4	BA-V-M4	BA-2-M4	BA-N-M4	BA-W-M4	BA-C-M4	BA-IA-M4
		DA-15O-M4-DN	DA-15A-M4-DN	DA-15V-M4-DN	DA-152-M4-DN	DA-15N-M4-DN	DA-15W-M4-DN	DA-15C-M4-DN	DA-15I-M4-DN
		DA-15UO-M4-DH	DA-15UA-M4-DH	DA-15UV-M4-DH	DA-15U2-M4-DH	DA-15UN-M4-DH	DA-15UW-M4-DH	DA-15UC-M4-DH	DA-15UI-M4-DH
1/8" NPT FEMALE		BA-O-F2	BA-A-F2	BA-V-F2	BA-2-F2	BA-N-F2	BA-W-F2	BA-C-F2	BA-IA-F2
		DA-15O-F2-DN	DA-15A-F2-DN	DA-15V-F2-DN	DA-152-F2-DN	DA-15N-F2-DN	DA-15W-F2-DN	DA-15C-F2-DN	DA-15I-F2-DN
		DA-15UO-F2-DH	DA-15UA-F2-DH	DA-15UV-F2-DH	DA-15U2-F2-DH	DA-15UN-F2-DH	DA-15UW-F2-DH	DA-15UC-F2-DH	DA-15UI-F2-DH

For Demand Check options add 'DC' to the end of the part number. For ISO colors replace the "U" with "I".

NPT Fittings

	FEMALE TEE	1/8" Female	NF-FT-F2
		1/4" Female	NF-FT-F4
	PIPE TEE	1/8" Female × 1/8" Female × 1/8" Male	NF-PTF2-F2M2
		1/4" Female × 1/4" Female × 1/4" Male	NF-PTF4-F4M4
	STREET TEE	1/8" Female × 1/8" Male × 1/8" Female	NF-STF2-M2F2
		1/8" Female × 1/8" Male × 1/8" Male	NF-STF2-M2M2
		1/4" Female × 1/4" Male × 1/4" Female	NF-STF4-M4F4
	90° ELBOW	1/8" Female × 1/8" Female	NF-EL90-F2F2
		1/4" Female × 1/4" Female	NF-EL90-F4F4
	90° STREET ELBOW	1/8" Female × 1/8" Male	NF-EL90-F2M2
		1/4" Female × 1/4" Male	NF-EL90-F4M4
	45° STREET ELBOW	1/8" Female × 1/8" Male	NF-EL45-F2M2
		1/4" Female × 1/4" Male	NF-EL45-F4M4
	FEMALE CROSS	1/8" Female	NF-CR-F2
		1/4" Female	NF-CR-F4
	HEX NIPPLE	1/8" Male × 1/8" Male	NF-HNM2-M2
		1/8" Male × 1/4" Male	NF-HNM2-M4
		1/4" Male × 1/4" Male	NF-HNM4-M4
	PIPE NIPPLE	1/8" Male, 1.5" Long	NF-PNM2-15
		1/8" Male, 2" Long	NF-PNM2-20
		1/8" Male, 2.5" Long	NF-PNM2-25
		1/8" Male, 3.5" Long	NF-PNM2-35
		1/4" Male, 1.5" Long	NF-PNM4-15
		1/4" Male, 2" Long	NF-PNM4-20
		1/4" Male, 2.5" Long	NF-PNM4-25
	COUPLER	1/8" Female × 1/8" Female	NF-COF2-F2
		1/4" Female × 1/8" Female	NF-COF4-F2
		1/4" Female × 1/4" Female	NF-COF4-F4
	CONNECTOR	1/8" Female × 1/8" Male	NF-CNF2-M2
		1/4" Female × 1/8" Male	NF-CNF4-M2
		1/4" Female × 1/4" Male	NF-CNF4-M4
	CAP	1/8" Female	NF-C-F2
		1/4" Female	NF-C-F4
	PLUG	1/8" Male	NF-P-M2
		1/4" Male	NF-P-M4
	HOSE BARB MALE	5/16" Hose Barb × 1/4" Male	NF-HB5-M4
		1/4" Hose Barb × 1/8" Male	NF-HB4-M2
		1/4" Hose Barb × 1/4" Male	NF-HB4-M4
	HOSE BARB FEMALE	1/4" Hose Barb × 1/8" Female	NF-HB4-F2
		1/4" Hose Barb × 1/4" Female	NF-HB4-F4
		5/16" Hose Barb × 1/4" Female	NF-HB5-F4



1111 Lakeside Drive
Gurnee, IL 60031-4099
+1 (866) 549-6446
www.ohiomedical.com

© 2020 Ohio Medical LLC. All rights reserved. This document contains information that is proprietary and confidential to Ohio Medical LLC. Use of this information is under license from Ohio Medical LLC. Any use other than that authorized by Ohio Medical LLC is prohibited. Ohio Medical and the Ohio Medical logo are registered trademarks of Ohio Medical LLC. Amvex, the Amvex logo, and Aeros are registered trademarks and Squire-Cogswell is a trademark of Ohio Medical LLC. Chemetron, Medstar and Oxequip are registered trademarks of Allied Healthcare Products, Inc. Ohmeda is a registered trademark of GE Healthcare Finland Oy Corporation. Schrader is a registered trademark of Schrader-Bridgeport International, Inc. Puritan Bennett is a registered trademark of Mallinckrodt Inc.