


Slide Bracket & Blank Wall Assembly



PART #	DESCRIPTION	MATERIALS OF CONSTRUCTION
261690	Wall Slide	Screws: Plated Steel
261394	Wall Cover Blank	Slide: Aluminum
261691	Console Slide	Cover Plate: Cycology®

Cycology® is a registered trademark of GE Plastics.

© 2018 Ohio Medical LLC.

This document contains information that is proprietary and confidential to Ohio Medical LLC. Use of this information is under license from Ohio Medical LLC. Any use other than that authorized by Ohio Medical LLC is prohibited. Ohio Medical LLC, Ohio Medical and the Ohio Medical LLC logo are registered trademarks of Ohio Medical LLC. Cycology is a registered trademark of Sabic Innovative Plastics IP B.V. Company.

Ohio Medical LLC | 1111 Lakeside Drive Gurnee, IL 60031 USA | 866.549.6446 | ohiomedical.com

255093 (Rev. 6) 04/2018