

## Cap and Chain Assembly



9/16" –18 Threads	Cap Size for:	Oxygen
3/4" –16 Threads	Cap Size for:	CO <sub>2</sub> , Medical Air, Nitrogen, N <sub>2</sub> O, & Vacuum
7/8" –14 Threads	Cap Size for:	Evacuation

<b>Part #</b>	<b>Description</b>
261683	Cap & Chain Assembly, CO <sub>2</sub> , Medical Air, Nitrogen, N <sub>2</sub> O, & Vacuum
261684	Cap & Chain Assembly, Evacuation
261685	Cap & Chain Assembly, Oxygen