

Outlet Surface Mounting Box

FEATURES & SPECIFICATIONS

- Single outlet surface mounting box (Latch valve and rough-in not included)
- Will accept any Ohio Medical® console outlet
- Four wall mounting holes $\frac{3}{16}$ " Diameter (bracket mounting screws provided by contractor)
- Dimensions: $3 \frac{1}{16}$ " W \times $2 \frac{7}{8}$ " D \times $4 \frac{7}{16}$ " H

PART NUMBERS

261636 Metal Cover Box & Metal Wall Bracket

261636-P Plastic Cover Box & Metal Wall Bracket



Latch Valve



Cover Box



Rough-in



Bracket with screws