



GAUGES

For Zone Valve Boxes and Isolation Valves

SPECIFICATION

Size	2"	
Style	ASME B40.1 Grade B	
Range	0 - 30" Hg	(0 - 700 mm Hg)
	0 - 100 PSI	(680 KPA)
	0 - 300 PSI	(2050 KPA)
Temperature	Ambient -13°F to 140°F (-25°C to 122°C)	
Media	Max 140°F	(122°C)
Connection	1/8" (6.35 mm) NPT Center Back	
Bourdon Tube	Copper Alloy	
Movement	Copper Alloy	
Dial	White Plastic with Stop Pin	
Pointer	Black Aluminum	
Case	Black Steel	
Window	Snap-in Acrylic	

Part Numbers

223016	High Pressure 0 - 300 PSIG
223017	Vacuum 0-30 Hg
223018	Low Pressure 0 - 100 PSIG



Vacuum



Low Pressure



High Pressure

© 2017 Ohio Medical LLC. All Rights Reserved.

This document contains information that is proprietary and confidential to Ohio Medical. Use of this information is under license from Ohio Medical. Any use other than that authorized by Ohio Medical is prohibited. Ohio Medical, and the Ohio Medical logo are registered trademarks of Ohio Medical LLC.

