





















Flowmeter Ordering Matrix - ISO

Instructions: Select the adapter, flowmeter and fitting to form an 11-digit part number.

Example: **7703-1280-921** **7703** DISS Nut & Gland **1280** Oxygen Flowmeter (0-15 LPM) **921** - Tubing Nipple

ADAPTERS	FLOWMETERS	FITTINGS
<p>No Adapter 7700 1/8" NPT Female</p>	<p style="text-align: center;">Oxygen Flowmeters</p> <hr/> <p>O₂ Flowmeter with O₂ DISS Fitting</p>  <p>1280 (0-15 LPM) (60 PSI) 1284 (0-15 LPM) 1109 (0 - 3.5 LPM) (60 PSI) 1319 (0 - 3.5 LPM) 1315 (0-1 LPM) 1318 (0 - 200 CCM)</p> <hr/> <p>O₂ Twin Flowmeter with O₂ DISS fitting</p>  <p>1292 (0-15 LPM)* (60 PSI)</p> <hr/> <p>O₂ Lo-High Duplex</p>  <p>1501 (LoFlo 3.5 & 0-15 LPM)* (60 PSI)</p> <hr/> <p style="text-align: center;">Medical Air Flowmeters</p> <hr/> <p>Air Flowmeter with O₂ DISS Fitting</p>  <p>1286 (0-15 LPM) (60 PSI) 1290 (0-15 LPM)</p> <hr/> <p style="text-align: center;">** "Y" adapter subject to change</p>	<p>No Adapter 907 DISS Fitting</p>
 <p>7701 Diamond</p>		 <p>921 Tubing Nipple Nut & Gland</p>
 <p>7702 DISS Hand-I-Twist (HIT)</p>	 <p>931 Clear Tubing Nipple Nut & Gland</p>	
 <p>7703 DISS Nut & Gland</p>		
 <p>7704 Schrader®</p>		
<p><i>Air Only</i></p>  <p>7705 Chemetron® Rectangle Striker (NCG)</p>		
 <p>7706 Chemetron® Round Striker (NCG)</p>		
 <p>7707 Puritan Bennett®</p>		
 <p>7708 O.E.S. Twist</p>		
 <p>7709 O.E.S. Med Star®</p>		
 <p>7714 DISS Male</p>		
 <p>7731 BOC Type</p>		
 <p>7732 DIN Type</p>		
 <p>7733 AFNOR Type</p>		
 <p>7734 AGA Type</p>		

© 2017 Ohio Medical, LLC. This document contains information that is proprietary and confidential to Ohio Medical, LLC. Use of this information is under license from Ohio Medical, LLC. Any use other than that authorized by Ohio Medical, LLC is prohibited. Ohio Medical and the Ohio Medical logo are registered trademarks of Ohio Medical, LLC. Chemetron and Med Star are registered trademarks of Allied Healthcare Products. Puritan Bennett is a registered trademark of Covidien LP. Schrader is a registered trademark of Schrader-Bridgeport International.

Flowmeters with Mounting Brackets for Dovetail Rail and Multi-Purpose Therapy Stand

Flowmeter and Manifold Bracket



Flowmeter with Hand-I-Twist (HIT)
DISS adapter, tubing nipple; manifold
bracket with one DISS male and two
DISS male checks
7721-1284-921 0-15 LPM O₂
7721-1290-921 0-15 LPM Air
7721-1319-921 LoFlo 3.5 O₂

Flowmeter and Manifold Bracket with Hose

Flowmeter with Hand-I-Twist (HIT)
DISS adapter, tubing nipple; manifold
bracket with one DISS male and two
DISS male checks, with 8' DISS HIT
female X DISS HIT female O₂ Hose
(not shown)
7700-0148-900

Universal Manifold Bracket

With four 1/8" NPT ports and
two plugs

0217-5359-800

Manifold Bracket With DISS Adapters

With two DISS male check units; one
DISS male adapter, one 1/8" NPT with
plug

6700-0411-800 O₂
6700-0412-800 Air



Tubing Nipples

Color-Coded Tubing Nipple

Constructed of Durable Plastic



Oxygen Wing Nut (O₂ DISS female)
& Gland with Tubing Nipple (*white*)
6700-0093-701 (pkg. of 20)
6700-0089-701 (1 ea)

Air Wing Nut (O₂ DISS male) &
Gland with Tubing Nipple (*black*)
6700-0093-702 (pkg. of 20)
6700-0089-702 (1 ea)

"Second Look" Clear Tubing Nipple
6700-0093-704 (pkg. of 20)
6700-0089-704 (1 ea)



DISS Union Wing Nut & Gland,
without chain, O₂ (*white*)
6700-0103-800



DISS Male with DISS Female N&G
(plastic)
6700-0725-800



Accessories



**8.5" Flowmeter or Vacuum
Regulator Security Chain**
with snap ring (stainless steel)
For securing equipment to
wall outlets.
6700-0404-200



**8.5" Flowmeter or Vacuum
Regulator Security Wire** with
secure hook. For securing
equipment to wall outlets.
6700-0403-200



O₂ Instant Flow Valve Fork
Flowmeter is "ON" but is turned
"OFF" by the bag on the fork.
• Immediate bag inflation
• Saves oxygen
6700-0400-200



O₂ Flow Selector Valve
Toggle switch selects flowmeter
O₂ or Humidified O₂.
6700-0402-200



Multi-purpose Rail Mount Slide
6700-0405-200



**Gauge Pole Clamp for .75"
- 1.5" Poles**, with 1/4" NPT
female and 1/8" NPT Tubing &
Nipple. Gas or Vacuum (shown)
6700-0406-200



Cylinder Wrench with Chain
6700-0020-300