

5700 Series Intermittent Suction Unit (ISO Product)

DESCRIPTION AND APPLICATION

The Ohio Medical® Intermittent Suction Unit is a combination continuous and intermittent regulator for use with hospital pipeline vacuum systems. Product has been an industry standard for decades, launched in 1975. It provides the versatility of an intermittent mode for nasogastric drainage and a continuous mode for tracheal, pharyngeal, or surgical suctioning. In the intermittent mode, the unit uses a pneumatic logic device to alternately turn the suction “ON” and “OFF.” This unique Unilogic Module sets the industry standard for providing a quiet and reliable means of cycling the regulator. The unit is equipped with independent controls to adjust the duration of the “ON” and “OFF” cycles, the level of vacuum, and the rate of flow.

The regulator also provides continuous vacuum, adjustable over the full range of available vacuum, to meet all clinical continuous suction requirements.

The Intermittent Suction Unit is classified as a high vacuum/high flow regulator.

FEATURES AND BENEFITS

- **Unilogic Module**

- Pneumatic

- Quiet, reliable operation

- Returns to atmospheric pressure during “OFF” cycle

- **Fully Adjustable**

- Timing cycles

- 3 to 30 seconds

- “ON” and “OFF” independently adjustable

- Factory preset at 15 seconds on/8 seconds off (begins in the off mode)

- **Flow rates**

- (Intermittent) 0 to 13 L/min, preset at 8 L/min (ASTM)

- (Continuous) 0 to 80 L/min (not adjustable)

- **Large, Easy to Read Gauge**

- Dual scale (kPa and mm Hg)

- 6.5 cm (2.5 in) gauge face

- Concave lens cover reduces glare

- Accuracy \pm 5% full scale

- **Modular Components**

- Simplifies service procedures

- Low cost replacement parts and kits



Intermittent Suction Unit (ISU)

- 0-200mm Hg
- 0-26 kPa

- **Durable Design and Construction**

- Few mechanical parts

- Rugged, shatter-resistant case

- **CE marked**

- **Full Three year parts and labor warranty**

- **Ordering Options**

Available in various configurations to include Probes for connection to source vacuum. Also available with options for collection devices and overflow safety protection.

SPECIFICATIONS*

Flow Rate

Intermittent – 0-13 L/min (preset at 8 L/min)
 Continuous – 0-80 L/min (not adjustable) without fittings at full increase

Timing

3 to 30 seconds
 Independent on and off
 Preset 15 seconds on/8 seconds off
 Begins in the "OFF" mode

Vacuum Range

0 to full range of available vacuum

Gauge Range

0-26 kPa (0-200 mmHg) and full vacuum

Gauge Accuracy

± 5% Full Scale Deflection

Dimensions

(H x W x D): 16.8 x 8.9 x 12.2 cm
 (6.6 x 3.5 x 4.8 in)

Weight

0.781 Kg (1 lb 11oz)

Warranty

Three Year

Standards

Conforms to ISO 10079-3 and ASTM F960

* Note: Specifications are nominal, subject to change without notice.

PURCHASING SPECIFICATIONS

Intermittent Suction Unit (ISU) with mode selection for intermittent, off and continuous suction. Intermittent mode utilizes a pneumatic logic device to cycle the unit automatically, returning the system to atmosphere in the "OFF" cycle and from 0 to Full Vac in the "ON" cycle. Timing cycles are independently adjustable from 3 to 30 seconds. Intermittent flow rate is adjustable from 0 to 13 Lpm, preset at 8 Lpm. Both mode selector and vacuum regulator controls are lockable. The regulator module and gauge assembly are interchangeable with the Continuous Vacuum Regulator.

- Dual Modes – intermittent and continuous
- Quiet, reliable operation
- Dependable, time tested gauge and controls
- Provides Full Vac regulation
- Low maintenance, easy to service
- Low cost replacement parts and kits

SHIPPING INFORMATION

Shipping weight:

1.1 Kg (2 lbs 7 oz) less fittings and adapters

Package Size:

(H x W x D): 20.3 x 8.9 x 12 cm
 (8 x 3.5 x 4.7 in)

Packaging:

One per carton

PART CONFIGURATOR

57XX - 1253 - 9XX

Intermittent Suction Unit

Adapter (vacuum source connection)

No Adapter with 1/8" NPT Female	00
British Std (BOC Type)	31
German Std *DIN Type)	32
French Std (AFNOR Type)	33
Nordic Std (AGA Type)	34

Fittings (patient connection)

No Fitting 1/8" NPT Female	00
Tubbing Nipple	01
Locking Gland	04
Locking Gland & Overflow Safety Trap	05
Locking Gland, Overflow Safety Trap and 1300 ml	06
Plastic Bottle with Wall Slide, Tubing	

© 2017 Ohio Medical, LLC.

This document contains information that is proprietary and confidential to Ohio Medical, LLC. Use of this information is under license from Ohio Medical, LLC. Any use other than that authorized by Ohio Medical, LLC is prohibited.

Ohio Medical and the Ohio Medical logo are registered trademarks of Ohio Medical, LLC.



Ohio Medical Authorized Representative

EC REP
 OxygenCare Ltd.
 2 Holfeld Business Park
 Kilmacanogue Co Wicklow Ireland
 Tel: +353 1 2769700 - Fax: +353 1 2764970

Ohio Medical, LLC
 1111 Lakeside Drive - Gurnee, IL 60031-4099 USA
 Toll Free: 866-549-6446 - Fax: 847-855-6218
 www.ohiomedical.com
 255499 (R2, 02/2017)