

instavac[®] C

MOBILE ASPIRATOR

Description & Application

Instavac C is a high vacuum & flow suction unit. It produces constant suction with vacuum & flow surpassing standard wall suction. Instavac C is packaged in a compact, durable molded Polystyrene case. The heart of the Instavac C is a fan cooled carbon rotary vane pump. Instavac C is insulated to provide one of the quietest sounding suction pumps on the market today. It is an excellent choice for the O.R., surgicenters, and doctors offices.

Features & Benefits

- Constant suction providing 550 mmHg and over 38 L/min flow
- Quiet operation
- Easy to clean and process for quick turnarounds
- Fan cooled pump, mechanical safety overflow trap and exhaust port filter

Ideal For

- G.I. Labs, Operating Rooms and Surgicenters requiring high vacuum and flow rates with quiet performance
- Surgical, scoping, and most specialty procedures



Part Number
757010

Shipping Information

Shipping Weight	48 lbs
Package Size	30" H × 24" W × 24" D
Packaging	One per carton



WITH RESPECT TO
ELECTRIC SHOCK, FIRE
AND MECHANICAL AND
OTHER HAZARDS ONLY IN
ACCORDANCE WITH UL
60601.1 AND CAN/CSA
C22.2 NO. 601.1

Technical Specifications, Instavac C 757010

Performance

Vacuum Range 0 to 550 mmHg

Free Air Flow 38 L/min minimum

Pump Fan cooled rotary carbon vane type

Collection Device

Canister Disposable plastic with mechanical shutoff

Capacity 1200cc standard

Tubing 14" with hydrophilic bacterial filter

Controls

Vacuum Regulator Rotary type on panel

Vacuum Gauge Calibrated 0 to 500 mmHg

Electrical Requirements

AC 120 volts, 60 Hz, 3A

Physical

Weight 37 lbs

Dimensions 28" H x 14" W x 16" D
(46" H with back pole)

Accessories



AI4120010
Disposable Collection
Canister, 1200cc (case of 10)



AI4120025
Disposable Collection
Canister, 1200cc (case of 25)



AI4125000
Ring Bracket for
1200cc Canister



AI5576C
Filter, Disposable Exhaust
(package of 25)



749002
Yankauer Suction Tip, Sterile



AI5785-14C
Filter, In-line suction with 14"
Tubing (case of 50)



749003
Suction Tubing with
End Connections,
72" x 1/4" ID, Sterile



756150
Dual Vacuum Port
Assembly Kit

© 2021 Ohio Medical LLC. This document contains information that is proprietary and confidential to Ohio Medical LLC. Use of this information is under license from Ohio Medical LLC. Any use other than that authorized by Ohio Medical LLC is prohibited. Ohio Medical and the Ohio Medical logo are registered trademarks of Ohio Medical LLC.