

# instavac® I

MOBILE ASPIRATOR

## DESCRIPTION AND APPLICATION

Instavac® I is a multipurpose portable suction unit. It produces both constant and intermittent suction with vacuum and flow levels equaling that of wall suction. Instavac® I is packaged in a compact, durable, molded Polystyrene case. The heart of the Instavac® I is a fan cooled carbon rotary vane pump. Instavac® I is insulated to provide one of the quietest sounding suction pumps on the market today.

## FEATURES AND BENEFITS

- Multipurpose unit, constant and intermittent suction modes
- Allows for standardization; replaces 3 single function style units
- Quiet operation
- Easy to clean and process for quick turnarounds
- Fan cooled pump
- Mechanical safety overflow trap and exhaust port filter
- 3 year warranty

## IDEAL FOR

Hospitals and extended care facilities requiring wall vacuum performance for Nasogastric, Thoracic, Oral, and Tracheal suctioning procedures.



## MULTI-PURPOSE SUCTION SYSTEM

### Constant or Intermittent Procedures

## TECHNICAL SPECIFICATIONS

### Performance

Vacuum Range: up to 380mmHg  
Free Air FLOW: 36 LPM minimum (constant mode)  
8 LPM (intermittent mode)  
Pump: fan cooled rotary carbon vane type

time Cycle (intermittent mode only)

Vacuum On: 20 seconds  
Vacuum Off: 10 seconds

### Collection Device

Canister: disposable plastic with mechanical shutoff  
Capacity: 1200cc standard  
Tubing: 14" with hydrophilic bacterial filter

### Controls

Vacuum Regulator: rotary type on panel  
Vacuum Gauge: calibrated to 300 mm Hg

### Electrical Requirements

AC: 120 volts, 60 Hz, 3A

### Physical

Weight: 34 lbs  
Dimensions: 28"H x 14"W x 16"D  
46"H with back pole

## Purchasing Specifications

Instavac®I Part Number: 756000

## Shipping Information

Shipping weight:  
48 lbs.

Package Size:  
30" H x 24" W x 24" D

Packaging:  
One per carton



WITH RESPECT TO ELECTRIC SHOCK,  
FIRE AND MECHANICAL AND OTHER  
HAZARDS ONLY IN ACCORDANCE WITH  
UL 2601-1 AND CAN/CSA  
C22.2 NO. 601-1 #9N56

## AVAILABLE ACCESSORIES

PART #	DESCRIPTION
756390	Disposable Chest Drainage Unit, Basket
756150	Dual Vacuum Port Assembly
AI4120010	Disp. Collection Canister, 1200cc (case of 10)
AI4120025	Disp. Collection Canister, 1200cc (case of 25)
AI2025	Bottle Cap Bracket for Slide Mount
AI2026	Bottle Cap Bracket for Aspirator Mount
AI2300C	Cap & Float Assembly, Disposable ( package of 25)
AI2311C	Turn Around Kit, Cap & Float, Tubing Filter (set of 10)
AI576C	Filter, Disposable Exhaust ( package of 25)
AI5785-14C	Filter, In-line suction w/14" Tubing (case of 50)
AI5785-18C	Filter, In-line suction w/18" Tubing (case of 50)



2017 Ohio Medical LLC. This document contains information that is proprietary and confidential to Ohio Medical LLC. Use of this information is under license from Ohio Medical LLC. Any use other than that authorized by Ohio Medical LLC is prohibited. Ohio Medical and the Ohio Medical logo are registered trademarks of Ohio Medical LLC.

