

# Hydrophobic Filter With Tubing Assembly

## **INSTRUCTIONS FOR USE**

1. Connect short tubing to vacuum outlet on unit.
2. Connect long tubing to collection canister. Use an elbow connector for a 800cc canister.
3. Replace filter per patient use or more often as needed.

**Reorder Number: 751101**



Ohio Medical Corporation  
1111 Lakeside Drive • Gurnee, Illinois • 60031-4099  
800-662-5822 • Fax 847-855-6218  
[www.ohiomedical.com](http://www.ohiomedical.com)

550118 (Rev.3) 01/07