

## tote-l-vac®

### DESCRIPTION AND APPLICATION

Tote-L-Vac® is a battery powered high vacuum and high flow suction unit. It features a soft, yet durable carrying case with easily accessible controls and battery. Although this unit is very compact, it provides high vacuum and high flow rates necessary for emergency oral aspiration procedures. Tote-L-Vac® accepts both 800cc and 1200cc standard collection canisters. Tote-L-Vac® is the perfect battery powered portable suction unit for EMS, home health, hospital patient transport and where ever battery power is a requirement.

### FEATURES AND BENEFITS

- Vacuum regulator and gauge
- Housed in a durable carrying case
- Easy to clean and process for quick turnarounds
- Compact and lightweight
- 3 year warranty

### IDEAL FOR

It is an excellent choice for EMS, hospital patient transport, extended care facilities, doctors offices, home healthcare or where ever a battery powered suction unit is required.



### BATTERY POWERED ASPIRATOR

E.M.S. • Patient Transport • Home Health Care • Surgicenters •  
Nursing Homes • Physician's offices

## TECHNICAL SPECIFICATIONS

### Performance

Vacuum Range: up to 550mmHg  
Free Air Flow: 36 LPM  
Pump: 12 VDC oil-less diaphragm type

### Collection Device

Canister: disposable plastic with mechanical shutoff  
Capacity: 800cc standard  
Tubing: 6' disposable sterile plastic  
Suction Tip: disposable Yankauer with control vent

### Controls

Vacuum Regulator: rotary type on panel  
Vacuum Gauge: Calibrated to 0-30"/760 mmHg

### Electrical Requirements

DC: 12 volts (nominal)

### Battery:

Type: rechargeable sealed lead-acid  
Capacity: 12 volts, 3.0 ampere-hour  
Average Run Time: 30 minutes when fully charged  
Charger: fast/float charger  
UL/1310 listed  
Charge Time: 8 hours or less to 85% charging,  
depending on initial charge of  
battery  
Charge Method: 110 volt outlet, or 12 VDC vehicle  
outlet with optional cord

### Physical

Weight: 15.5 lbs  
Dimensions: 8"H x 12"W x 9.5"D

## AVAILABLE ACCESSORIES

PART #	DESCRIPTION
AI5100125	12VDC Cigarette Power Cord
749146	Replacement Lead Acid Battery
749121	AC Charger
AI4080010	Disp. Collection Canister, 800cc (case of 10)
AI4080025	Disp. Collection Canister, 800cc (case of 25)
750150	Filter, Hydrophobic with 5" Tubing (case of 50)

## PURCHASING SPECIFICATIONS

tote-I-vac® Part Number: 4900000

## SHIPPING INFORMATION

Shipping weight:  
20 lbs.

Package Size:  
12" H x 8" W x 8" D

Packaging:  
One per carton

Terms:  
net 30 days

ISO 13485



CERTIFICATE  
NO. FM 39489



TRUSTED BRANDS OF OHIO MEDICAL™

aeros®  
Division of Ohio Medical

AMVEX Corporation  
Division of Ohio Medical

Ohio  
Medical Corporation®

Squire-Cogswell  
Division of Ohio Medical

Ohio Medical Corporation  
1111 Lakeside Drive - Gurnee, IL 60031-4099 USA  
TollFree: 866-549-6446 - Fax: 847-855-6218  
www.ohiomedical.com  
Form No. 550150 (Rev.10) 10/2011