

Introducing: Sure-Trap™

Overflow Safety Trap

The new, cost-effective **Sure-Trap** is designed to save you time and money, while also protecting your medical suction equipment. With the new **Sure-Trap**, there is no cleaning required--simply discard when the trap is contaminated. The **Sure-Trap** is the last line of defense and protection against accidental contamination of your vacuum regulators, wall outlets and central pipeline systems.



Features

- ◆ Float design ensures maximum vacuum flow and positive shutoff
- ◆ Leak-Free cap is permanently bonded to the jar
- ◆ Lightweight materials are latex free
- ◆ 1/8" NPT fitting for direct attachment to vacuum regulator
- ◆ Helps healthcare professionals comply with:
 - NFPA 99: Health Care Facilities 2015
 - DOD UFC 4.510.01.2016

Sure-Trap with 1/8"x1/8" NPT Resin Fitting	
Package	Order Number
12 traps	7725-0365-912
25 traps	7725-0365-925
50 traps	7725-0365-950

To place an order, contact your Sales Representative, or call 1.866.549.6446