

Sure-Trap™

Overflow Safety Trap

The new, cost-effective Sure-Trap™ is designed to save you time and money, while also protecting your medical suction equipment. With the new Sure-Trap™, there is no cleaning required—simply discard when the trap is contaminated. The Sure-Trap™ is the last line of defense and protection against accidental contamination of your vacuum regulators, wall outlets and central pipeline systems.



FEATURES

- Float design ensures maximum vacuum flow and positive shutoff
- Leak-Free cap is permanently bonded to the jar
- Lightweight materials are latex free
- 1/8" NPT fitting for direct attachment to vacuum regulator
- Helps healthcare professionals comply with:
 - NFPA 99: Health Care Facilities 2015
 - DOD UFC 4.510.01.2016

Sure-Trap with 1/8"x1/8" NPT Resin Fitting

Package	Order Number
12 traps	7725-0365-912
25 traps	7725-0365-925
50 traps	7725-0365-950

To place an order, [contact your Sales Representative](#), or call 1.866.549.6446



EMERGO EUROPE
Prinsessegracht 20 2514 AP The Hague The Netherlands

AUSTRALIAN SPONSOR: EMERGO AUSTRALIA
Level 20, Tower II Darling Park 201 Sussex Street Sydney, NSW 2000 Australia

