

## Amvex Disposable Hydrophobic Suction Filter

### DESCRIPTION & APPLICATION

The Amvex Hydrophobic Filter offers protection for Vacuum Regulators and Vacuum Systems from contamination. These easy to use in-line filters have a 0.2 µm PTFE membrane designed to protect against an over-flow of the collection canister and aqueous aerosols. The Amvex Disposable Hydrophobic Filters are constructed to using a durable Borosilicate and PTFE membrane in a Polypropylene Housing.

### FEATURES & BENEFITS

- Retains aqueous aerosols preventing Vacuum Regulator repairs
- Retains bacteria and viruses
- Suitable for use with portable suction equipment to protect pump
- Bacteria filtration efficiency 99.98%

#### Tubing Nipple Version



1/4" Hose Barb x 1/4" Hose Barb

VR-FLTR-HB4HB4

#### Threaded Version



1/4" Hose Barb x 1/8" Male NPT

VR-FLTR-M2HB4