

Gas Regulator & Flowmeter with Humidifier



Pin Index Cylinder Regulator Option
GR-MU-00101



Bullnose Cylinder Regulator Option
GR-MU-00100

Gas Regulator and Flowmeter with Humidifier

Product Benefits



Flowmeter

- Provides accurate gas flow measurement and control within a range of 1-15 liters per minute
- For use in a variety of respiratory therapy clinical applications
- Features large easy to read increments with an expanded 1-5 L/min range for improved readability at low flows
- Dual taper and rear surface textured to enhance readability
- Supply pressure of 414 kPa/60 psi
- Flush setting delivers a minimum of 65 L/min
- Lifetime warranty on flowtube & shroud



Gas Regulator

- Customizable to meet your facility's specific needs
- 100% tested and quality inspected
- Calibrated to 60 PSI
- Machined brass



Humidifier

- Lightweight and compact
- Shatter-resistant jar and lid
- Chrome-plated brass connections
- Quiet, reliable operation
- Service parts available
- Holds 300 mL
- Incorporates pressure relief valve