



# Gas Regulators

## Click Style and Compact Gas Regulators


### Product Benefits

- ▶ Fully customizable to meet your facility's specific needs
- ▶ Choose regulator style, connection and gas type
- ▶ ISO or USA color options
- ▶ Click style regulator has a built in flowmeter
- ▶ 100% tested and quality inspected




Visit [www.amvex.com](http://www.amvex.com) for the gas regulator part configurator and more!

# Gas Regulators Click and Compact style



## Click Style

- ① Amvex® Click Style is available in 50 PSI Preset
- ② Available with a CGA Nut and Nipple, CGA Yoke Style & Bullnose
- ③ Available in Oxygen and Air
- ④ Aluminum body with brass core
- ⑤ Built in Flowmeter available in 0-4 LPM, 0-8 LPM, 0-15 LPM, and 0-25 LPM increments



## Compact Style

- ① Amvex Compact Regulators are available in 50 or 60 PSI Preset
- ② Available with a CGA Nut and Nipple, CGA Yoke Style and Bullnose
- ③ Available in Oxygen, Air, N<sub>2</sub>O, CO<sub>2</sub>, Nitrogen and Heliox
- ④ Available in USA and ISO colors
- ⑤ Available with standard Thorpe Style Flowmeter options
- ⑥ Machined brass

(\*MR conditional option available)

## Part Configurator

Basic matrix shown. For customization or power take off options, please contact your Amvex Sales Representative. Not all configurations are possible. All configurations are 50 PSI Preset unless specified.

### Click Style Regulators

## GR(X)X-YYYGX-X

<b>Gas Regulator Type:</b> <b>Flow Range (litres per minute)</b> 4 LPM: <b>4</b> Mini 15 LPM: <b>M1</b> 8 LPM: <b>8</b> Mini 25 LPM: <b>M2</b> 15 LPM: <b>15</b> Mini 8 LPM: <b>M8</b> 25 LPM: <b>25</b> Mini 4 LPM: <b>M4</b>	<b>Connection Type:</b> YOKE: <b>870</b> NUT & NIPPLE: <b>540</b> BULL NOSE OXY: <b>BNO</b>	<b>Optional Colors:</b> Green: <b>G</b>	<b>Flow Port:</b> Nipple: <b>T</b> DISS Male: <b>D</b>
			<b>Label:</b> Amvex: <b>X</b>

### Compact Regulators

## GR-WWW-(XXXX)-(YY)(Z)

*Compact Gas Regulators*  
 NOTE: For ISO option, please add "I" after the part number

<b>Gas and Connection Type:</b> Oxygen Nut & Nipple: <b>540</b> N <sub>2</sub> O Nut & Nipple: <b>326</b> Oxygen Yoke Style: <b>870</b> N <sub>2</sub> O Yoke Style: <b>910</b> Medical Air Nut & Nipple: <b>346</b> CO <sub>2</sub> Nut & Nipple: <b>320</b> Medical Air Yoke Style: <b>950</b> CO <sub>2</sub> Yoke Style: <b>940</b> NIT Yoke Style: <b>960</b> NIT Nut & Nipple: <b>580</b> CGA 540 1/4 FNPT OXY: <b>F40</b> Heliox Nut & Nipple: <b>280</b> CGA 346 1/4 FNPT AIR: <b>F4A</b> OXY Preset Bull nose: <b>BNO</b>	<b>Optional Flowmeter:</b> 8 LPM Flowmeter: <b>FM08</b> 16 LPM Heliox: <b>FM16</b> 12 LPM: <b>FM12</b> 15 LPM Dial: <b>FMD1</b> 15 LPM Flowmeter: <b>FM15</b> 25 LPM Dial: <b>FMD2</b> 15 LPM 60PSI: <b>FMN1</b>	<b>Additional Options:</b> ISO: <b>I</b> ISO/MRI: <b>L</b> ISO/B/KPA 50 PSI: <b>N</b> ISO/60 PSI: <b>H</b>
		<b>Additional Outlet:</b> Power Take-off For FM: <b>PT</b> Nothing: <b>XX</b> Swivel for FM: <b>SW</b>

© 2016 Ohio Medical Corporation. This document contains information that is proprietary and confidential to Ohio Medical Corporation. Use of this information is under license from Ohio Medical Corporation. Any use other than that authorized by Ohio Medical Corporation is prohibited. Ohio Medical Corporation and the Ohio Medical logo are registered trademarks of Ohio Medical Corporation. Amvex and the Amvex logo are registered trademarks of Ohio Medical Corporation.

TRUSTED BRANDS OF OHIO MEDICAL®



Amvex®  
 25B East Pearce Street, Richmond Hill, ON., L4B 2M9 Canada  
 Toll Free: 866-462-6839 - Fax: 905-764-7743

www.amvex.com

MR-GR-SPEC-I 01/16