

# Amvex® Gas Regulators


## CLICK STYLE

**Built-in Flowmeter**  
available in 0-4 L/min, 0-8 L/min,  
0-15 L/min, and 0-25 L/min increments



## COMPACT STYLE

Available with  
standard Thorpe Style Flowmeter options



### Click Style Regulators

## GR(X)X-YYYGX-X

Gas Regulator Type		Connection Type		Flow Port	
4 L/min	4	Yoke	870	Nipple	T
8 L/min	8	Nut & Nipple	540	DISS Male	D
15 L/min	15	Bullnose O <sub>2</sub>	BNO		
25 L/min	25				

- Customizable to meet your facility's specific needs
- Choose regulator style, connection and gas type
- Click style regulator has a built in flowmeter
- 100% tested and quality inspected

### Compact Regulators

## GR-WWW-(XXXX)-(YY)(Z)

Gas & Connection Type		Optional Flowmeter	
O <sub>2</sub> Nut & Nipple	540	8 L/min	FM08
O <sub>2</sub> Yoke Style	870	12 L/min (N <sub>2</sub> O & CO <sub>2</sub> )	FM12
Medical Air Nut & Nipple	346	15 L/min	FM15
Medical Air Yoke Style	950	15 L/min (60 psi)	FMN1
N <sub>2</sub> Nut & Nipple	580	16 L/min (Heliox)	FM16
N <sub>2</sub> Yoke Style	960	15 L/min (Dial)	FMD1
N <sub>2</sub> O Nut & Nipple	326	25 L/min (Dial)	FMD2
N <sub>2</sub> O Yoke Style	910		
CO <sub>2</sub> Nut & Nipple	320		
CO <sub>2</sub> Yoke Style	940		
Heliox Nut & Nipple	280		
O <sub>2</sub> CGA 540 1\4" FNPT	F40		
Air CGA 346 1\4" FNPT	F4A		
O <sub>2</sub> Preset Bullnose	BNO		
Air Preset Bullnose	BNA		

Additional Outlet	
Power Take-off for FM	PT
Swivel for FM	SW

Additional Options	
ISO	I
ISO (60psi)	H

Basic matrix shown. For customization or other options, please contact your Ohio Medical Sales Representative. Not all configurations are possible. All configurations are 50 psi preset unless specified.