

big ben[®] Blood Pressure Monitor

Used with ri-former wall diagnostic stations 97-LF3655-109.018

Product Specifications

- Glass-fiber reinforced wall bracket for maximum strength (comparable to steel). Simple conversion to the rail model
- Generously dimensioned scale Ø 5.8" or 5.7" and high-contrast scale layout for optimum readability
- Spacious basket on the reverse side for tubes and cuff
- Spiral tube can be extended to 8'
- Supplied with latex-free Adult Cuff (1)
- Non-wearing air-release valve with fine adjustment
- Specially hardened copper-beryllium membrane, resistant to aging, can withstand pressures of up to maximum 600 mmHg
- Micro-filter to protect the release valve and measurement system
- No zero-point fixation required
- Maximum error tolerance: ± 2 mmHg
- BHS validation, see: bihsoc.org/bp-monitors/

Additional Cuff Sizes Available


97-LF108.004	NIPB Cuff, Big Ben, Adult
97-LF121.002	NIPB Cuff, Big Ben, Adult SM
97-LF123.004	NIPB Cuff, Big Ben, Adult XL
97-LF128.001	NIPB Cuff, Big Ben, Thigh XL
97-LF130.003	NIPB Cuff, Big Ben, Child
97-LF133.002	NIPB Cuff, Big Ben, Infant
97-LF135.002	NIPB Cuff, Big Ben, Newborn
97-LF126.002	NIPB Cuff, Big Ben, Thigh

WALL MODEL 97-LF1465.018

Can be swiveled horizontally 65° to the left and right



MOBILE MODEL 97-LF1478.018



- Height adjustable from 35.4" to 55.1"
- Diameter 20.9"
- Sturdy mobile base made of fiberglass and chrome plated pole
- Easy to move due to low friction castors, anti-skid wheels, 2 locking brakes
- Spacious basket on the reverse side (included)



Optional pole basket (97-10392)