

April 24, 2024

Venturi Backpack End of Life

Dear Valued Customer,

This letter is to inform you that effective immediately, the Ohio Medical Venturi Backpack Suction Accessory has been discontinued. Affected part numbers are listed below.

Discontinued Part Numbers for Venturi Backpack: 6700-0690-800, 6702-0690-907



Venturi Backpack

See next page for discontinued part numbers of suction regulators with venturi backpack.

As an alternative, we would like to recommend our battery operated portable suction unit that we believe will satisfy your requirements. Our Care-e-vac[®] 3 (758000) is a portable continuous, battery-operated suction device that may also be combined with a wheeled stand (AI2580000 and adapter kit 758007) for enhanced mobility.

While we regret any inconvenience this announcement may cause, we are eager to support your future needs with our quality product lines. If you have any questions, please reach out to our customer service team or your local Ohio Medical sales representative.

Regards,

Suzanne Moyer

Suzanne Moyer
Sr. Business Development Manager
Ohio Medical



Suction Regulators with Venturi Backpack

Ohio Brand Suction Regulator Part Numbers

Product numbers with 2, 3, 4 or 5 in indicated position below.

XX XX - XXXX - 9 XX

REGULATOR TYPE

CVR-ISU 67

PTS CVR-ISU 87

PUSH-TO-SET VACUUM REGULATORS

	Analog		Digital	
	Air	Oxygen	Air	Oxygen
Three Mode CVR	2224	3224	4224	5224
Two Mode CVR	2225	3225	4225	5225
Three Mode High CVR	2226	3226	4226	5226
Two Mode Low CVR	2227	3227	4227	5227
Two Mode High CVR	2330	3330	4330	5330
Two Mode Neo CVR	2231	3231	4231	5231
Two Mode Ped CVR	2234	3234	4234	5234
ISU	2251	3251	4251	5251
Ped ISU	2271	3271	4271	5271
Three Mode High CVR	2229	3229	4229	5229
Two Mode High CVR	2230	3230	4230	5230
Two Mode Neo CVR	2232	3232	4232	5232
Two Mode Ped CVR	2233	3233	4233	5233
ISU	2253	3253	4253	5253
Ped ISU	2273	3273	4273	5273

PATIENT CONNECTION

No Fitting with 1/8" NPT Female 00

Tubing Nipple 01

Locking Gland & Overflow Safety Trap 05

DISS Male & Overflow Safety Trap with DISS Wing Nut & Gland 08

WALL CONNECTION

No Adapter with 1/8" NPT Female	00
Ohmeda	01
DISS Handtight	02
Chemetron (Air)	05
Chemetron (Oxygen)	06
Puritan Bennett	07
British BOC	31
German DIN	32
French AFNOR	33
Scandinavian AGA	34

Amvex Brand Suction Regulator Part Numbers

VR-SSTU-VWX-(Y)(Z) where Z = F, U or V