



SCORE Quick Guide

Surgery Center Operating Room Equipment



From Starting Up

Whether you're building from the ground up or retrofitting a current facility, you want reliable products so you won't lose precious time when you're caring for patients. Ohio Medical carries a variety of medical gas and suction and oxygen therapy products to meet your needs.

To start, Ohio Medical has a vast array of medical gas equipment to fit your space and keep your operating costs to a minimum.

Vacuum Systems

The oil-sealed duplex vacuum systems and oil-less rotary vane vacuum systems come in a variety of sizes, operate quietly and come with touch screen panels. They're also compact, which is important when space is at a premium.

Manifolds & Alarms

To keep track of how well the systems are working, we also provide manifolds and master and area alarms. Their touch screen panels make them easy to operate and read in a variety of lighting conditions.

Outlets


In addition, we also provide a variety of outlets to meet your facility's specifications. They are offered in wall, ceiling and console installation. Standard hose drops with retractors and custom length hoses are readily available.



Rotary Vane Vacuum Systems



Rotary Claw Vacuum Pumps




Area Alarms



Manifolds



Hose Drops



Outlets

To Gearing Up

We also carry a full line of suction and oxygen therapy equipment when you're ready to gear up all your clinical pre-op, post-op or operating room needs.

Vacuum Regulators

We carry a complete line of suction regulators and accessories designed to provide reliable clinical solutions, ease-of-use and patient safety.

Integrated Flowmeters & Flowmeters

Our flowmeters meet strict standards of durability and precision. And our patented integrated flowmeter has one unique outlet and two reliable ports.

Diagnostic Equipment

Offering includes thermometers, blood pressure devices and more to help you measure and observe various aspects of your patient's health. Available in wall mount and portable options.

Portable Suction

We also have portable suction solutions in the event of a power outage. The battery-operated units are compact and ready to use when you need them.

Hose Assemblies & Rail Systems

Customize your hose and rail systems for your facility. Ohio Medical has the ability to develop these products to your specifications.


When you're ready to build or replace your equipment, our dedicated team of professionals will take care of you to ensure you receive the right products for your facility. Call 866-549-6446 today to get your quote.



Vacuum Regulators



Diagnostic Equipment



Portable Suction



Rail Systems



Integrated Flowmeters & Flowmeter



Hoses

To Fixing Up

If your products need servicing, go directly to the source to get them fixed! We offer a variety of service and maintenance programs for your Ohio Medical products that will meet your budget and scheduling requirements.

Medical Gas Equipment Repairs

Our Preventative Maintenance Contracts for Vacuum, Air and Medical Gas include:

- Oil and filter changes
- Inspection and adjustment of belts, drive couplings, mountings, etc.
- Test performance, measure amp draws
- Recommend necessary repairs, service parts, and upgrades
- Gas Manifolds: check for leaks, verify proper operation, line pressure and adjust regulators
- Gas Alarms: test audible alarms, verify visible alarms and pressure settings

Service contract customers receive a 10% discount on oil, filters, and repair parts. Labor for repairs outside the standard PM's during the term of the contract will also receive a 10% discount from our standard rate.

To schedule Service call 847-855-6251 and for Technical Support call 847-855-6234

Medical Device Equipment Repairs

We also provide maintenance for your Medical Devices. We're opening a new service center dedicated to taking care of all your repair needs from fixing a cracked faceplate on a vacuum regulator to replacing a knob on a flowmeter. And, if we can't fix it, we manufacture the real thing, so replacements are easy!

Items we repair include:

- Vacuum Regulators
- Portable Suction
- Integrated Flowmeters
- Flowmeters



To submit a repair request visit, [ohiomedical.com/contact-us/returns](https://www.ohiomedical.com/contact-us/returns)