

# Integrated Flowmeter


Thorpe Tube Style



# ONE UNIQUE OUTLET | TWO RELIABLE PORTS

For New Construction and Retrofit Applications

Ohio Medical has a flowmeter that won't be moved or lost. The Patented Tube Style Integrated Flowmeter combines a tube style flowmeter and a medical gas outlet into one simple, compact design – helping to reduce replacement/repair costs. The gas outlet offers an extra connection – giving you the flexibility to attach another flowmeter, hose or ventilator. Now you can do more, in less space. Rest assured you will always have a flowmeter when you and your patients need one most.



## Respiratory Therapists & Nurses:

*Reduce Worry & Frustration*

Stop searching for a flowmeter. With the Integrated Flowmeter, once it's installed, it will always be there when you need it.

## Biomedical Engineers & Department Managers:

*Avoid Unnecessary Replacement/Repair Costs*

The Tube Style Integrated Flowmeter can help you save time and money because it reduces the number of lost or damaged flowmeters you reorder and/or repair each year.

## Equipment Planners & Design Engineers:

*Save Space*

The Tube Style Integrated Flowmeter can be customized to your health care institution's needs and saves space with its overall compact design.

## Facilities & Maintenance:

*Reduce Maintenance Costs*

Stop moving the flowmeter from outlet to outlet, room to room, wearing out the outlets and causing them to leak. Because the Tube Style Integrated Flowmeter is permanently affixed to the wall, it stays put and lowers maintenance costs.

## AUXILIARY PORT OPTIONS



DISS




Chemetron



Diamond/Ohmeda

The Patented Integrated Flowmeter comes equipped in DISS, Chemetron and Diamond/Ohmeda connections giving you the flexibility to deliver more solutions.

# SPECIFICATIONS



## IMPROVE SPEED & SAFETY

Once installed, the Integrated Flowmeter is permanently affixed to the wall. No more looking for flowmeters. This flowmeter also comes with a 50 psi outlet to connect to other devices.

## ELIMINATE REPAIRS OR REPLACEMENTS

This wall-secured flowmeter reduces outlet damage due to frequent removal of traditional flowmeters. Outlet and flowmeter replacement costs are lowered because the Integrated Flowmeters cannot be removed.



0-200 cc/min and 0-1 L/min models carry a 1 year warranty on the flowmeter portion

# AUXILIARY PORT APPLICATIONS



Use the extra connection to add another flowmeter or connect a hose to another device.

# INTEGRATED FLOWMETER CONFIGURATOR FOR RETROFIT APPLICATIONS

Available in different flowrates, gas types and styles of connectors – the Integrated Flowmeter is retrofit capable and fits many different outlet types. Refer to the Integrated Flowmeter Ordering Guide (FI-OGUIDE) for additional information.

## FI-TXXYZ-WN-S

FLOWRATE	
200 cc/min	<b>2C</b>
1 L/min	<b>01</b>
3.5 L/min	<b>03</b>
8 L/min	<b>08</b>
15 L/min	<b>15</b>

COLOR CODE	
USA	<b>U</b>
ISO	<b>I</b>

GAS TYPE	
Oxygen	<b>O</b>
Medical Air	<b>A</b>

AUXILIARY PORT	
DISS Male	<b>D</b>
Diamond/Ohmeda	<b>O</b>
Chemetron	<b>C</b>

OUTLET COMPATIBILITY	
Diamond 1	<b>A*</b>
Diamond 2 & 3 (Gemini)	<b>B†</b>
Hill-Rom BeaconMedaes DiamondCare (Extended Barrel Version)	<b>C</b>
Ohio Medical Amico BeaconMedaes Series B Gentec	<b>E</b>
Chemetron Series 500, 560 (NCG)	<b>F*</b>
Chemetron Series 148, 400, 460 (NCG)	<b>G*</b>
Tri-Tech (Ohmeda Style)	<b>T*</b>

\*Requires a Retrofit Kit | †For Wall Mount Only

# INTEGRATED FLOWMETER OUTLETS FOR NEW CONSTRUCTION

## 261XXX-Y(Y)-TZZ

AUXILIARY PORT	
DISS	<b>000</b>
Diamond/Ohmeda	<b>010</b>
Chemetron	<b>020</b>

GAS/OUTLET TYPE	
Oxygen/Wall	<b>1</b>
Oxygen/Console	<b>2</b>
Air/Wall	<b>13</b>
Air/Console	<b>14</b>

FLOWRATE	
200 cc/min	<b>2C</b>
1 L/min	<b>13</b>
3.5 L/min	<b>03</b>
8 L/min	<b>08</b>
15 L/min	<b>15</b>



+



+



= 261010-1-T15

For ISO colors add '-ISO'  
(Example: 261010-1-T15-ISO)

Back Body/Rough-In

Trim Plate

Integrated Flowmeter

© 2019 Ohio Medical, LLC. This document contains information that is proprietary and confidential to Ohio Medical, LLC. Use of this information is under license from Ohio Medical, LLC. Any use other than that authorized by Ohio Medical, LLC is prohibited. The Ohio Medical logo is a registered trademark of Ohio Medical, LLC. Chemetron, MedStar and Oxequip are registered trademarks of Allied Healthcare Products, Inc. Ohmeda is a registered trademark of GE Healthcare Finland Oy Corporation. Puritan Bennett is a registered trademark of Mallinckrodt Inc. Schrader is a registered trademark of Schrader-Bridgeport International, Inc. Hill-Rom is a registered trademark of Hill-Rom Services Inc.